

recoveriX

NEUROTECHNOLOGY


www.recoveriX.at

Micah Ching
G.tec Neurotechnology Hong Kong Limited



Brain-Computer Interfaces for Stroke Rehabilitation

recoveriX



Companies

- Guger Technologies OG, Graz, Austria
- g.tec medical engineering GmbH, Schiedlberg, Austria
- g.tec medical engineering Spain SL, Barcelona, Spain
- g.tec neurotechnology USA Inc., Albany, USA
- g.tec neurotechnology Hong Kong Inc., HK





Research Projects

- **H2020 Eurostars: seizureAI** – epilepsy monitoring and deep learning with EEG
- **H2020 Eurostars: ComaAid** – combination of EEG and fNIRS for DOC patients
- **H2020 Eurostars: EEG-DDS** – combination of EEG with decision support system
- **FFG: BrainGait** – combination of BCIs for rehabilitation with treadmills
- **H2020: MultiSense** – combination of DBS with physiology
- **FFG: recoveriX-Leg** – lower limb stroke rehabilitation with BCI
- **H2020: Rhumbo** – magnetic stimulation of the cortex



Research Projects

- H2020: HOPE – high frequency oscillation detection in epilepsy
- H2020 ITN: Pro-Gait – combination of BCIs with exoskeletons
- H2020 ITN: DOCMA – international training network for BCIs and DOC
- Marie-Curie: MoveAgain – combination of recoveriX with TMS, tDCS
- H2020: Astonish – fNIRS and EEG
- H2020 Eurostars: ComAlert – coma prediction
- EC Flagship: Graphene – development of Graphene electrodes


Application Examples

MEASURE BIODIVERSITY RELIABLY - EVEN IN THE HIMALAYAS

During an Austrian expedition to Chulu Far West in Nepal (Camp III), g.NEObits was used to measure the effects of high altitude on both EEG and ECG parameters. The expedition started in Dhaulagiri at an altitude of 7000 m near Annapurna I. The team gained between 300 and 600 m each day and settled in base camp at 4800 m. After one night in base camp, the high camp was established at 5600 m on the Chulu glacier. At 3 AM, the team started to climb Chulu Far West and reached the summit at 11 AM. g.NEObits was used to record EEG channels over sensor motor areas and ECG channel from 2 expedition members. The persons performed a self-paced finger movement every 10 seconds. The onset and offset of each movement was recorded by an external switch connected to g.NEObits.

INTERACT WITH ROBOTS

The Serbian artist Dragan Fic equipped a KUKA robot with hundreds of pencils, and used it to create numerous works of art with BCI software and g.NEObits. He selects drawing commands just by thinking, and the robot paints on a vertical and a horizontal wall panel by Dragan's mind. The performance was shown during the AIS Electronics Festival in Ljubljana, and visitors could also control the device.



Application Examples

The first Cybathlon took place on October 8th, 2016 at the ETH Zurich. People with disabilities competed side by side at the Brain-Computer-Interface Race (BCI Race). Pilot Numa Pospjoly from the Team BrainTweakers from Ecole Polytechnique Fédérale de Lausanne (EPFL) in Switzerland won Gold in the BCI Race using his brain waves to control an avatar along a race-track on a virtual train. The game is called "BrainRunners" and was specially developed for the Cybathlons BCI Race. The Brain Tweakers have been using high-quality g.tec equipment.

g.USBamp Research Edition



g.USBamp RESEARCH is the non-certified version of g.USBamp. Therefore, it is less expensive and intended to be used for research applications only. The device has freely accessible and configurable features, although any kind of ECoG experiments are excluded.

User Experience

Natalie Mrachacz-Kersting, PhD, Associate Professor
Aalborg University, Denmark

"I have used the g.USBamps for approximately four years now, and the quality of the signals surpasses our greatest expectations. Last year, we also invested in the g.Neotilus system, and have been able to move from more stringent laboratory conditions to real world settings, where BCIs are becoming more prominent. I would highly recommend both systems to anyone involved in BCI research."



13




Recommended setup for a fully equipped BCI lab



g.tec medical engineering GmbH, office@gtec.at, www.gtec.at, tel. +43 7251 22240, fax +43 7251 22240 39

Check products of interest and request an offer per fax (+43 7251 22240 39) or e-mail (office@gtec.at)!

Brain-Computer-Interface (BCI)



"A system for controlling a device e.g. computer, wheelchair or a neuroprosthesis by human intention which does not depend on the brain's normal output pathways of peripheral nerves and muscles" [Wolpaw et al., 2002].


- HCI – Human Computer Interface
- DBI – Direct Brain Interface (University of Michigan)
- TTD – Thought Translation Device (University of Tübingen)

recoveriX




Physiological background

Imagination of hand movement causes an ERD which is used to classify the side of movement. The desynchronization occurs in motor and related areas of the brain.



recoveriX



Stroke Rehabilitation

- Motor imagery (MI) based rehabilitation was proven to be an effective therapy.
- BCI based stroke rehabilitation combines a physical therapy with BCI-based motor imagery decoding into an integrative rehabilitation strategy




How recoveriX works

We combine:

- Motor imagery
 - Functional electrical stimulation
 - Visual feedback via the movement of virtual hands (Mirror therapy)
- Simultaneous stimuli (Activity of the motor cortex + actual movement of the hand + visual feedback)
 - Plasticity of the brain is supported to recover motor functions
 - Reduce negative symptoms

recoveriX




How recoveriX works

- Measures EEG signals and instructs the patient to imagine a left/right hand movement
- BCI detects if the corresponding motor cortex is activated
 - Trigger a electric stimulation
 - Display movement of an avatar hand (VR) = Mirror therapy
- Movement of the hand + Mirror neurons activates the sensorimotor cortex
 - More activation in the brain

recoveriX





How recoveriX works



recoveriX
NEUROTECHNOLOGY




recoveriX setup




recoveriX setup


EEG gel to the electrodes



FES electrodes to Extensor digitorum communis muscle





Start setting the stimulation with the healthy side!



recoveriX


Videos

EEG montage, FES placement, therapy – [video](#)




recoveriX

FES electrode placement




FES electrodes to Extensor digitorum communis muscle

Start setting the stimulation with the healthy side!



recoveriX


rX paradigm





Seconds

0
1
2
3 Sound
4 Cue: Right or Left
5
6 Feedback
7
8

Visual (Avatar)
FES



recoveriX



Read more about recoveriX

- recoveriX Website: www.recoveriX.at
- recoveriX Results: www.recoveriX.at
- recoveriX Patient Interviews: www.youtube.at/gtecmecaleengineering

25
Sep
2020
recoveriX results: huge functional improvement in chronic stroke patient
Schwabing Upper Austria | Read more

19
Nov
2020
recoveriX results: Motor improvement 30 years after stroke
Innsbruck | Read more

18
Sep
2020
recoveriX results: Stroke patient improves dramatically after no success in the first 17 sessions
Innsbruck | Read more

recoveriX

